

**PRIMAL CONTACTS - Soft (Hydrophilic) Contact Lenses**  
**Eye Care Instructions**  
**IMPORTANT: Please read carefully and keep this information for future use.**

Primal Brand Contacts Lenses are available as corrective and non-corrective lens designs. The Lenses are composed of 62% Polymacon and 38% water a silicone-containing hydrogel.

**PRIMAL (NON-CORRECTIVE) CONTACT LENSES PARAMETERS**

Diameter : 14,0 mm to 14,8 mm  
Base Curve : 8,6 mm to 8,8 mm  
Power : Plano

**INDICATIONS (USES) AND WEARING RESTRICTIONS**

Primal® Brand Fashion Lenses is indicated for non-corrective use in the colour enhancement of the iris and have no therapeutic effects for daily wear.

**CONTRAINDICATIONS (REASONS NOT TO USE PRIMAL BRAND CONTACT LENSES)**

Inflammation or infection in or around the eye or eyelids.  
Any eye disease, injury, or abnormality that affects the cornea, conjunctiva, or eyelids.  
Any previously diagnosed condition that makes contact lens wear uncomfortable.  
Any systemic disease that may affect the eye or be exaggerated by wearing contact lenses.  
Allergic reactions of ocular surfaces of the eye (adnexa) that may be induced or exaggerated by wearing contact lenses or use of contact lens solutions.  
Allergy to any ingredient, such as mercury or Thimerosal, in a solution which is to be used to care for your contact lenses.  
Severe dry eye.  
Any active corneal infection (bacterial, fungal, protozoal or viral).  
Not to be used if you are unable to follow the lens care regimen or unable to obtain assistance to do so.

**POTENTIAL ADVERSE SIDE EFFECTS FROM WEARING PRIMAL BRAND CONTACT LENSES**

Stinging, burning, itching (irritation), pain, tearing or other eye discomfort or sensation of foreign object in the eye.  
Excessive watering, unusual eye secretions, redness of the eye or difficulty opening after sleep.  
Poor vision, blurred vision, rainbows or halos around objects, sensitivity to light (photophobia).

Should side effects occur due to the use of **Primal®** Brand contact lenses, remove the contact lens immediately, and consult your Eye Care Professional immediately and before resuming in the use of contact lenses.

**WARNINGS**

- **Problems with contact lenses or lens care products could result in serious injury to the eye.** Proper use and care of your contact lenses and lens care products, including lens cases, are essential for the safe use of these products.
- Eye problems, including a sore or lesion on the cornea (corneal ulcers) can develop rapidly and lead to loss of vision.
- The risk of an infected sore or lesion on the cornea (ulcerative keratitis) is greater for individuals who wear contact lenses past the recommended time in the instructions or from your Eye Care Professional. **(Do not wear lenses for longer than 8 hours a day)**
- When Individuals wear their lenses overnight (outside the approved use), the risk of an infected sore or lesion on the cornea (ulcerative keratitis) is greater than among those who do not wear them while sleeping.
- If you experience eye discomfort, excessive tearing, vision changes, redness of the eye or other problems, you should immediately remove your lenses and promptly contact your Eye Care Professional.
- It is recommended that you see your Eye Care Professional routinely as directed.

**STORAGE AND VALIDITY**

- Storage: Seal and store the lenses in the proper lens container with saline soaking solution and keep at room temperature (15~20°).
- Validity: Unopened - 5 years from the date of manufacturing. Opened – 3 months from the time of opening.

## APPLYING AND REMOVING YOUR CONTACTS

### APPLYING STEP BY STEP :

**Thoroughly wash your hands with soap and dry using a clean towel.**

**Place the lens on the tip of your right index finger. Be sure the lens is positioned correctly.**

**Use your other hand to gently hold your upper eyelid from blinking. Pull down your lower eyelid with the other fingers of your applying hand.**

**Look up at the ceiling and gently place the lens on the lower part of your eye.**

**Slowly release your eyelid and close your eye for a moment. Then blink a few times to center the lens.**

### REMOVING STEP BY STEP

**Thoroughly wash your hands with soap and dry using a clean towel.**

**Look up at the ceiling and gently pull down holding your lower eyelid.**

**Using your index finger gently slide the lens to the lower white part of your eye. Slowly release your eyelid and close your eye for a moment. Then blink a few times to center the lens.**

**Gently squeeze the lens between your thumb and index finger and remove it.**

### CHEMICAL LENS DISINFECTION

Clean the contact lenses and case with a recommended cleaning solution every 2 days (even if not used).

Thoroughly rinse lenses with a recommended rinsing solution.

Store lenses in a thoroughly clean contact lens case in a recommended storage solution.

Leave lenses in unopened storage case until use. Rinse thoroughly with fresh solution prior to next wear.

### CARE REGIMEN AND SAFETY PRECAUTIONS

#### HANDLING PRECAUTIONS

- **DO NOT** use if the sterile blister package is opened or damaged.
- **CAREFULLY** follow the handling, insertion, removal and wearing instructions provided with your contact lenses and those prescribed by your Eye Care Professional.
- **ALWAYS** wash and rinse hands thoroughly before handling lenses. Do not get soaps, lotions, cosmetics, sprays or deodorants in contact with your eyes or lenses.
- **ALWAYS** apply makeup after contacts have been put in your eyes.
- **ALWAYS** handle lenses gently and avoid dropping them.
- **DO NOT** touch the lens with finger nails.
- **NEVER** use tweezers or other tools to remove your lenses from the lens container unless specially indicated for that use. Pour the lens and solution directly in to the palm of your hand.

#### LENS WEARING PRECAUTIONS

- **IMMEDIATELY** remove your lenses if they become red or irritated or if any other side effects as described before occur.
- **DO NOT DRIVE** while wearing Primal Brand costume lenses. They may obscure vision.
- **NEVER** wear contact lenses beyond the amount of time recommended by the Eye Care Professional. Contact lenses should not be worn for longer than 8 hours a day as this will increase the risk of infection or permanent damage.

**DO NOT** sleep while wearing lenses. It is recommended to remove lenses during a period of rest overnight or longer.

Ask your Eye Care Professional about wearing contact lenses during sports activities.

**AVOID** going in the water while wearing lenses. This includes lakes, the ocean, swimming pools, hot tubs and even the shower. The water may contain unseen threats to your eyes.

**AVOID** environments that are extremely dry, dusty, handling of harmful chemical substances, exposure to moisture, or extreme changes in temperatures. This may cause irritation or infections of the eye. Consult your Eye Care Professional if you are subject to these environments.

**AVOID** using contact lenses if you suffer any eye diseases or are allergic to materials contained in the lens.

**ALWAYS** inform your employer of being a contact lens wearer. Certain jobs require the use of eye protection or may require the individual not wear contact lenses.

## **SOLUTION PRECAUTIONS**

**ALWAYS** use the recommended product for cleaning, rinsing and storage of contact lenses. Different lens care products cannot always be used together.

**ALWAYS** use fresh, unexpired lens care solutions.

If the lens sticks (stops moving or cannot be removed) apply 2 to 3 drops of the recommended lubricating or wetting solution directly to the eye and wait until the lens begins to move freely before it is removed.

**USE ONLY** a chemical (not heat) lens care system. Use of heat (thermal) care systems can damage your contact lens.

**DO NOT** use saliva or anything other than the recommended solutions for lubricating or wetting lenses.

**ALWAYS** make sure you completely cover your lenses with the recommended storage solution when not wearing the lenses (stored). If your lenses start to dry out it will make them more difficult to become wet again.

## **LENS CASE PRECAUTIONS**

Contact lens cases should be emptied, cleaned, and rinsed with recommended solution and allowed to air dry. Lens cases should be replaced at regular intervals as recommended by your Eye Care Professional.

### **Topics to discuss with Your Eye Care Professional :**

For proper use of **Primal**® Brand soft contact lenses and eye care maintenance consult your Eye Care Professional.

**FOLLOW UP VISITS** with you Eye Care Professional are necessary to assume the continuing health of your eyes.

**AVOID** using contact lenses if you suffer any eye diseases or are allergic to materials contained in the lens.

**ALWAYS** contact your Eye Care Professional before using any medicine in your eyes.

## **THEATRICAL AND COSMETIC CONTACT LENSE WAIVER AND RELEASE OF LIABILITY**

**Cosmetic contact lenses are like other contact lenses, but are intended solely to change the normal appearance of the eye in a decorative/cosmetic fashion.**

**These products are intended to be worn only with the assistance of a Qualified Eye Care Professional and professional eye care supervision. Do not wear these lenses if you have not had an eye exam in the past year, have never worn contact lenses or have not been trained in contact lens insertion, removal and lens care.**

**Like other contact lenses, cosmetic contacts can cause a variety of eye injuries and conditions. Contact lenses can cause infection, redness, corneal ulcer and/or other problems if not used properly. If left untreated, certain ocular (eye) problems such as infection, can lead to corneal scarring which can lead to vision impairment. In extreme cases, corneal ulcer can result in blindness and eye loss. Other risks include conjunctivitis, corneal edema, allergic reaction, abrasion, and reduction in visual acuity, contrast sensitivity and other visual functions resulting in interference with driving and other activities.**

**These risks cannot be sufficiently controlled unless:**

- 1. The wearer obtains contact lens wearing advice from an eye care professional.**
- 2. The lenses are evaluated by an eye care professional.**
- 3. The wearer remains under appropriate professional supervision, including undergoing yearly eye exams.**

**I certify that I have read this document and I fully understand its content. I am 18 years of age or older.**

**I agree to follow all of the instructions listed above which includes wear and care instructions for contact lenses.**

**I am aware that this is a waiver and release of liability and I sign it of my own free will.**

**I agree to not hold Streets Management Group Pty Ltd t/a Shindigs or its directors or its staff liable for any injury to my eyes or for any financial loss that may occur directly or indirectly from the use of cosmetic contact lenses.**

**Name:** \_\_\_\_\_ **(Please print clearly)**

**Signature:** \_\_\_\_\_ **Date:** \_\_\_\_\_

**ID: Supplied Form of ID:** \_\_\_\_\_ **ID Number :** \_\_\_\_\_